

Electrostatic Enhanced Compact Aerosol Concentrator Development

**Presented at the 2003 Joint Service Scientific
Conference on Chemical & Biological Defense
Research**

Peter Coyle*, Tim Pletcher & Chetna Bindra

Sarnoff Corporation

2003-11-19



Report Documentation Page				Form Approved OMB No. 0704-0188	
Public reporting burden for the collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to Washington Headquarters Services, Directorate for Information Operations and Reports, 1215 Jefferson Davis Highway, Suite 1204, Arlington VA 22202-4302. Respondents should be aware that notwithstanding any other provision of law, no person shall be subject to a penalty for failing to comply with a collection of information if it does not display a currently valid OMB control number.					
1. REPORT DATE 19 NOV 2003		2. REPORT TYPE N/A		3. DATES COVERED -	
4. TITLE AND SUBTITLE Electrostatic Enhanced Compact Aerosol Concentrator Development				5a. CONTRACT NUMBER	
				5b. GRANT NUMBER	
				5c. PROGRAM ELEMENT NUMBER	
6. AUTHOR(S)				5d. PROJECT NUMBER	
				5e. TASK NUMBER	
				5f. WORK UNIT NUMBER	
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES) Sarnoff Corporation				8. PERFORMING ORGANIZATION REPORT NUMBER	
9. SPONSORING/MONITORING AGENCY NAME(S) AND ADDRESS(ES)				10. SPONSOR/MONITOR'S ACRONYM(S)	
				11. SPONSOR/MONITOR'S REPORT NUMBER(S)	
12. DISTRIBUTION/AVAILABILITY STATEMENT Approved for public release, distribution unlimited					
13. SUPPLEMENTARY NOTES See also ADM001851, Proceedings of the 2003 Joint Service Scientific Conference on Chemical & Biological Defense Research, 17-20 November 2003. , The original document contains color images.					
14. ABSTRACT					
15. SUBJECT TERMS					
16. SECURITY CLASSIFICATION OF:			17. LIMITATION OF ABSTRACT UU	18. NUMBER OF PAGES 23	19a. NAME OF RESPONSIBLE PERSON
a. REPORT unclassified	b. ABSTRACT unclassified	c. THIS PAGE unclassified			


Electrostatic Collection Technology

Acknowledgements

- **The Electrostatic Enhanced Compact Aerosol Concentrator Development Program is part of the Technology Transfer Program being funded by DTRA / SBCCOM, under the technical direction of the Edgewood Chemical Biological Center**
- **Technical Point of contact Dr. J. Bottiger, Dr. E Steubing**

Electrostatic Collection Technology


...versatile technology capable of interfacing to many detection modalities...



Electrostatic Enhanced Compact Aerosol Concentrator Development


...Roadmap...

FY 03




**Experimental
Electrostatic Capture
Module**

FY 04



**Inline Aerosol
Concentrator**


FY 05



**Compact Electrostatic
Aerosol Concentrator
Prototype**



Hydrophobic Membrane



**Wetted Column Electrostatic
Capture Module**

Electrostatic Enhanced Compact Aerosol Concentrator Development


...Description of Current Effort...

- **Moderately concentrated aerosol drawn through experimental capture duct**
 - simulates the output of a cyclone concentrator
 - monodispersed aerosol consisting of fluorescent beads
 - sheath air employed to confine aerosol stream
- **Electrostatically charged by corona array**
- **Particles electrostatically focused directly into small liquid volume**
 - beads deposited on solid surface in initial experiments
- **Component configuration modified to optimize collection efficiency**
- **Performance determined by fluorescent imaging techniques**

Electrostatic Enhanced Compact Aerosol Concentrator Development


...experimental set-up for optimizing electrostatic capture configuration...

Orthogonal Collection Approach




Electrostatic Focusing of Particles Into a Liquid Using Corona Charging & Transport

Operation Schematic for Orthogonal Collection Approach



Electrostatic Enhanced Compact Aerosol Concentrator Development


...experimental set-up for optimizing electrostatic capture configuration...




Compact Electrostatic Aerosol Concentrator

...compact, low power, high flow rate...

- **Direct aerosol concentration without energy consumptive inertial separation process upstream**
 - Pressure drop orders of magnitude lower than inertial separation collector
 - 7 unit array samples 210 LPM with 1 watt fan
- **Particles charged by radial array of corona tips**
- **Particles deposited onto small liquid wetted column located on the axis of the aerosol duct**




Electrostatic Focusing of Particles Onto Liquid Wetted Column Concept



Electrostatic Focusing of Particles Onto Column Experiments

Test Setup




2.5 mm diameter
nozzle

Corona Array

1.5 mm dia. dry
column

ITO coated glass
focus electrode




Electrostatic Focusing of Particles Onto a Column Experiments

Electrostatic deposition of smoke onto a 1.5 mm diameter dry column @ 1 liter per minute aerosol flow rate through 2.5 mm diameter nozzle

ITO coated glass
focus electrode


1.5 mm diameter dry column
ground electrode




- Excess deposition beyond column occurs when deposited particles form insulating layer preventing newly arriving particles from discharging to the column ground
- The particles then repel each other causing them to be deposited on the surrounding surface.
- This will not occur in the actual implementation of the technology where the aerosol density is orders of magnitude lower than that of the smoke used to visualize the aerodynamic characteristics of the collector.

Electrostatic Focusing of Particles Onto Column Experiments – AEROSOL DATA


Electrostatic deposition of 2.3 μ beads onto a 2 mm diameter dry column with 3 mm diameter nozzle @ 1 liter per minute aerosol flow rate






Electrostatics Off



Electrostatics On



**Electrostatics Off 10X
Electrostatics On 10X**




> 99% of 2.3 μ Particles Removed From Aerosol by Electrostatic Collection


Electrostatic Focusing of Particles Onto Column Experiments – AEROSOL DATA

Electrostatic deposition of room air onto a 2 mm diameter dry column with 3 mm diameter nozzle @ 1 liter per minute aerosol flow rate



> 90% of 0.5 μ m Particles removed from Aerosol by Electrostatic collection



Electrostatics Off




Electrostatics On



Electrostatic Focusing of Particles Onto Column Experiments – DEPOSITION DATA


Electrostatic deposition of 2.3 μ fluorescent beads onto a 2 mm diameter dry column with 3 mm diameter nozzle @ 1 liter per minute aerosol flow rate




Top of Column

Deposition predominantly on column collection surfaces

Small amount on the focusing electrode surface



Side of Column



Focus Electrode

Half this density of particles observed on focusing electrode surface with the electrostatics turned off.

Corona Charging Technology

...developed on independent program involving electrostatic deposition of pharmaceuticals ...

■ Corona charging advancements

- Developed new corona charging technique superior to currently available technology
- High efficiency (demonstrated >99%) uni-polar charging
- Current controlled corona eliminates ion density variation due to tip erosion
- Array of Corona tips generates uniform corona

■ Corona wind management


- “wind” blows particles toward collection port
- controlled by tip array

Charging


- Corona-charging
- Optimum Charge Saturation

Corona Charging Technology

- **Advantages of corona charging**
 - Predictable particle charging
 - Charging obeys Pauthenier equation
 - Field charger
 - Uses maximum E-Field for capture
 - Designed for theoretical efficiency for respirable particle size



Corona Tip Array



Charging

- Corona-charging
- Optimum Charge Saturation

Self Cleaning Liquid-Air Interface


...withstand differential pressure with efficient particle transport...

■ Teflon treated nylon mesh retains liquid in collection channel

- Debris tolerant
- 40% open area
- Low cost - standard Teflon fabric treatment
- thermal imbedment assembly

■ Self Cleaning

- over pressurize channel to purge liquid past mesh to rinse contaminants off surface
- liquid seal re-established when pressure is reduced



*Results in clog-free,
low-maintenance
continuous operation*

Liquid

- Self-cleaning
Hydrophobic
Membrane

Electrostatic Surface Deposition

...Benefits ...

- **Superior to impaction methods**
 - **Better control of deposition parameters**
 - Smaller size sample
 - Potential to deposit mono layer
 - **Less likely to damage pathogen**
 - Gentle Low Velocity Deposition
 - **Lower pressure drop**
 - Lower power
 - Quieter aerosol collection fan
- **Tape transport from collector to detector**
 - **Rapid, automated, precise location**
 - **Tape configuration:**
 - Disposable - reel to reel
 - Reusable - loop
 - Materials - polymer, ceramic, metal

Surface

- Low Velocity Deposition

Electrostatic Sorting

...sample enrichment by electrostatic deflection of selected particles...

■ Description of operation

- Stream of particles electrostatically charged by corona created at the exit of nozzle
- Particles interrogated by UV laser
- resulting laser induced fluorescence used to determine particles of biological origin
- Selected particles focused onto small diameter target
- rejected particles deflected to annular surface surrounding target

■ Features



- High S/N ~ 50,000:1
- High switching speeds > 10Khz
- High efficiency (>99% @ 2um)

Surface


Sorting

Electrostatic Sorting

...feasibility experiment...



Experimental result showing deposition of HP cyan toner (3um-9um) particles onto collection electrode



Experimental result showing deposition of HP cyan toner (3um-9um) particles onto rejection electrode


Experimental electrostatic sorting mechanism

Surface

Sorting

Electrostatic Sorting

...deposition onto tape concept ...



Surface

Sorting

Benefits of Electrostatic Collection

- **Affords plug and play detector integration for standardized air sampling module**
 - Deposition into liquid for assay and antibody based detection
 - Solid surface deposition for spectrographic detection
 - Electrostatically focused aerosol stream for advanced particle sorting techniques
- **Electrostatically focus particles into μl liquid volumes**
 - Interfaces to advanced detection technologies
 - Clog free hydrophobic membrane
 - Capable of low temperature operation
- **Electrostatic surface deposition superior to impaction methods**
 - Better control of deposition parameters
 - Less likely to damage pathogen
 - Lower pressure drop reduces fan power & noise